


## VWR Collection compact incubators

### VWR Collection INCU-Line digital incubators IL 10 & IL 23



#### Compact design

- **Digital PID control with PT100 sensor and LED display**
- **Temperature range: ambient +5 °C to 70 °C**
- **Excellent temperature stability**
- **Easy to use, safe and reliable**

The VWR digital incubators IL 10 & IL 23 offer an economical and space saving solution for microbiology or haematology applications. The temperature can be freely adjusted in increments of 0,1 °C up to a maximum of 70 °C. Heating elements are incorporated into the housing base and walls to ensure temperature stability and uniformity. The temperature control can be re-adjusted with a certificated reference temperature measuring instrument.

The housing and door frame are constructed entirely of epoxy coated mild steel. The VWR IL 10 is available with a transparent plexiglass window to provide full visibility or a painted, opaque plexiglass window to protect photosensitive samples.

The VWR IL 23 is supplied with two perforated, stainless steel shelves which can be placed in three different positions. The VWR IL 10 is supplied with one perforated shelf and can, optionally, be loaded with up to two specially designed aluminium tube racks for vertical incubation of a maximum of 24 dip slide tubes.

Model	IL 10	IL 23
Interior volume	10 l	23 l
Temperature range	ambient +5 to 70 °C	
Temperature accuracy	± 0,5 °C	
Temperature stability (at 37 °C)	± 0,5 °C	
Exterior dimensions w x d x h (mm)	285 x 280 x 335	340 x 360 x 400
Interior dimensions w x d x h (mm)	230 x 220 x 200	290 x 290 x 270
Weight (Kg)	8,3	12,8
Voltage	230 V / 50 Hz	

Description	Cat. No.
Mini Incubator IL 10, digital, transparent window	390-0384
Mini Incubator IL 10, digital, painted opaque window	390-0385
Incubator IL 23, digital, transparent window	390-0482
Extra shelf for IL 10, stainless steel, 226 x 215 mm	390-0387
Extra shelf for IL 23, stainless steel, 286 x 285 mm	390-0483
Tube rack for dip slides, 2 x 6 positions, Ø 32 mm	390-0386

## Your European Distribution Partner

### Austria

VWR International GmbH  
Graumanngasse 7  
1150 Wien  
Tel.: 01 97 002 0  
Fax: 01 97 002 600  
E-mail: [info@at.vwr.com](mailto:info@at.vwr.com)

### Belgium

VWR International bvba  
Researchpark Haasrode 2020  
Geldenaaksebaan 464  
3001 Leuven  
Tel.: 016 385 011  
Fax: 016 385 385  
E-mail: [customerservice@be.vwr.com](mailto:customerservice@be.vwr.com)

### Denmark

VWR - Bie & Berntsen  
Transformervej 8  
2730 Herlev  
Tel.: 43 86 87 88  
Fax: 43 86 87 90  
E-mail: [info@dk.vwr.com](mailto:info@dk.vwr.com)

### Finland

VWR International Oy  
Valimotie 9  
00380 Helsinki  
Tel.: 09 80 45 51  
Fax: 09 80 45 52 00  
E-mail: [info@fi.vwr.com](mailto:info@fi.vwr.com)

### France

VWR International S.A.S.  
Le Périgares – Bâtiment B  
201, rue Carnot  
94126 Fontenay-sous-Bois cedex  
Tel.: 0 825 02 30 30 (0,15 EUR TTC/min)  
Fax: 0 825 02 30 35 (0,15 EUR TTC/min)  
E-mail: [info@fr.vwr.com](mailto:info@fr.vwr.com)

### Germany

VWR International GmbH  
Hilpertstrasse 20a  
D - 64295 Darmstadt  
Tel.: 0180 570 20 00\*  
Fax: 0180 570 22 22\*  
E-mail: [info@de.vwr.com](mailto:info@de.vwr.com)  
\*0,14 €/Min. aus d. dt. Festnetz,  
Mobilfunk max. 0,42 €/Min.

### Hungary

VWR International Kft.  
Simon László u. 4.  
4034 Debrecen  
Tel.: (52) 521-130  
Fax: (52) 470-069  
E-mail: [info@hu.vwr.com](mailto:info@hu.vwr.com)

### Ireland / Northern Ireland

VWR International Ltd /  
VWR International  
(Northern Ireland) Ltd  
Orion Business Campus  
Northwest Business Park  
Ballycoolin  
Dublin 15  
Tel.: 01 88 22 222  
Fax: 01 88 22 333  
E-mail: [sales@ie.vwr.com](mailto:sales@ie.vwr.com)

### Italy

VWR International s.r.l.  
Via Stephenson 94  
20157 Milano (MI)  
Tel.: 02 332 03 11  
Fax: 800 152 999  
E-mail: [info@it.vwr.com](mailto:info@it.vwr.com)

### The Netherlands

VWR International B.V.  
Postbus 8198  
1005 AD Amsterdam  
Tel.: 020 4808 400  
Fax: 020 4808 480  
E-mail: [info@nl.vwr.com](mailto:info@nl.vwr.com)

### Norway

VWR International AS  
Haavard Martinsens vei 30  
0978 Oslo  
Tel.: 02290  
Fax: 815 00 940  
E-mail: [info@no.vwr.com](mailto:info@no.vwr.com)

### Poland

Labart Sp. z o.o.  
A VWR International Company  
Limbowa 5  
80-175 Gdansk  
Tel.: 058 32 38 210  
Fax: 058 32 38 205  
E-mail: [labart@labart.pl](mailto:labart@labart.pl)

### Portugal

VWR International - Material de  
Laboratório, Lda  
Edifício Neopark  
Av. Tomás Ribeiro, 43- 3 D  
2790-221 Carnaxide  
Tel.: 21 3600 770  
Fax: 21 3600 798/9  
E-mail: [info@pt.vwr.com](mailto:info@pt.vwr.com)

### Spain

VWR International Eurolab S.L.  
C/ Tecnología 5-17  
A-7 Llinars Park  
08450 - Llinars del Vallès  
Barcelona  
Tel.: 902 222 897  
Fax: 902 430 657  
E-mail: [info@es.vwr.com](mailto:info@es.vwr.com)

### Sweden

VWR International AB  
Fagerstagatan 18a  
163 94 Stockholm  
Tel.: 08 621 34 00  
Fax: 08 621 34 66  
E-mail: [info@se.vwr.com](mailto:info@se.vwr.com)

### Switzerland

VWR International AG  
Lerzenstrasse 16/18  
8953 Dietikon  
Tel.: 044 745 13 13  
Fax: 044 745 13 10  
E-mail: [info@ch.vwr.com](mailto:info@ch.vwr.com)

### UK

VWR International Ltd  
Customer Service Centre  
Hunter Boulevard  
Magna Park  
Lutterworth  
Leicestershire  
LE17 4XN  
Tel.: 0800 22 33 44  
Fax: 01455 55 85 86  
E-mail: [uksales@uk.vwr.com](mailto:uksales@uk.vwr.com)

Go to [www.vwr.com](http://www.vwr.com) for the latest news, special offers and details of your VWR distribution partner